PoloDent LC Flowfill is a visible light cure easy flowing radiopaque microhybrid composite for anterior and limited posterior restorations and for sealing. It is based on methacrylate-resin and inorganic filler particles of 0.05-1 µm. The total filler load is 62% and the total filler volume 47%. Radiopacity of LC Flowfill is 150% aluminium (radiopacity of enamel: ~200% aluminium, of dentine: ~100% aluminium). The composite material meets the requirements of DIN EN ISO 4049, type 1, class 2, group 1.

Indications:

- Restorations of class III, IV and V
- Small class I and II restorations
- Lining of cavities
- Sealings

Side Effects:

Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection areas close to the pulp should be covered with a thin layer of pulp capping material.

Incompatibility with other materials:

Do not use in combination with substances containing eugenol because eugenol inhibits the polymerization of the composite. Neither store the composite material in proximity of eugenol containing products, nor let the composite allow coming into contact with materials containing eugenol.

Preparation of tooth and color matching:

Clean the tooth with a fluoride-free polishing paste (e.g. flour of pumice with water) prior to preparation and color matching.

Ascertain the tooth shade while teeth are still moist and select the appropriate LC Flowfill shade.

LC Flowfill is shaded according to VITA® -Shades.

Cavity preparation:

After isolation (best with a rubber dam) prepare the cavity with minimal tooth reduction.

For pulp protection areas close to the pulp should be covered with a thin layer of pulp capping material.

Condition and apply a bonding agent according to the manufacturer instructions. The total etch technique and subsequent application of a light cure bonding agent is recommended. Otherwise use a self etching bonding agent before application of LC Flowfill.

Placement of LC Flowfill:

Apply LC Flowfill in the selected shade. Place it by using a suitable instrument. Matrix strips may be used.

Application from syringe:

Screw off the cap of the syringe and substitute it by one of the supplied Application tips. After use discard the application tip and attach the cap to the syringe again. Application tips are for single use only due to hygienic reasons.

Curing of LC Flowfill:

For an optimum result apply LC Flowfill in layers of max. 1 mm. Light cure each increment separately.

By using a polymerization unit (wavelength range 400-500 nm) with a light intensity of at least 1000mW/cm², cure each increment as follows:
Lighter shades (e.g. A1, A2, B2) 20 seconds
Darker shades (e.g. A3.5, B3) 30 seconds

Hold the light emission window as close as possible to the filling material.

Finishing:

Depending on the shape of the area contour with finishing diamonds, flexible polishing discs or carbide burs. Polish to high gloss with polishing discs or polishing brushes.

Storage:

Do not store above 25°C (77°F). Protect from direct sunlight. Do not use after expiry date. Close syringe immediately after use to avoid exposure to light.

Applicator Tips are for single use only due to hygienic reasons.

Additional notes/warnings:

- The ambient light of the dental lamp may start polymerization of the composite.
- Do not use any resin to adjust viscosity of composite restorative material.
- Avoid contact with skin, mucous membrane and eyes.
- Unpolymerized composite may have an irritant effect and can lead to sensitization against methacrylates.
- Commercial medical gloves do not protect against the sensitizing effect of methacrylates.
- VITA® is a registered trade mark of the VITA-Zahnfabrik, Bad Säckingen, Germany.

Warranty:

DS Dental Supply GmbH warrants this product will be free from defects in material and manufacture. DS Dental Supply GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and DS Dental Supply GmbH's sole obligation shall be repair or replacement of the DS Dental Supply GmbH product.

Limitation of Liability:

Except where prohibited by law, DS Dental Supply GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

**Keep away from children!
For dental use only!**

Any serious incident related to this medical device must be reported to the manufacturer and the competent authorities.

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