Endo Stand

INSTRUCTIONS FOR TREATMENT

en	Endo Stand	03
nl	Endo standaard	04
fr	Support pour endo	05
de	Endo-Stand	06
pl	Stojak Endo	07
da	Endo stand	08
sv	Endo-stativ	09
no	Endo stativ	10
fi	Endo-teline	11

POLO | DENT | [®]

E N D O

Product #	Description
50 20 00	Endo Stand

 **Polo MB Oisterwijk BV**
Laarakkerweg 5
5061 JR Oisterwijk, The Netherlands

 **Larident S.r.l.**
Via Lamaneigra 12R
16030, Tribogna (GE) Italia

MANUFACTURER'S INSTRUCTIONS FOR THE TREATMENT OF MEDICAL DEVICES ACCORDING TO THE EN ISO 17664:2017 STANDARDS

WARNING NOTICE:

THE PRODUCT IS NOT SUPPLIED STERILE. STERILIZE IT BEFORE USE.

Autoclavable plastic support for endo sponges

Restriction on treatments: the product can be reused after autoclaving it

Area of use: support for endodontic instruments

Storage and transport: the product does not have specific needs

Method of use: the instrument is a plastic support for endodontic sponges. It's hollow inside and it can contain disinfectant liquids for providing support to the sponge during the cleaning of root canal endodontic instruments (files and reamers).

Disinfection and sterilization: sterilize in an autoclave machine before use at 134° for 4 minutes.

Maintenance: no specific requirement

Disposal: dispose of the product according to the correct disposal procedures

Function check and test: visual check for any damage or deformation of the product

Storage: the product does not have specific needs

The instructions listed above have been validated by the manufacturer of the medical devices as suitable for preparing the device for use.

It is the user's responsibility to ensure that the preparation process, the equipment the materials used and the staff obtain the desired results.

Any deviation from the recommended instructions must be carefully checked with regards to its efficiency and possible adverse consequences.

Any serious incident related to this medical device must be reported to the manufacturer and the competent authorities.

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