Scaling Tips

STERILIZATION INSTRUCTIONS

en	Scaling tips	03
nl	Scaling tips	05
fr	Embout de détartrage	07
de	Ultraschallspitze	09
pl	Końcówka do skalingu	11
da	Scaler spid	13
sv	Scalerspets	15
no	Skaleringstips	17
fi	Hammaskivenpoistokärki	19

 Polo MB Oisterwijk BV
Laarakkerweg 5
5061 JR Oisterwijk, The Netherlands

 Guilin Yikeshi Medical Instrument Co., Ltd.
D-8, Guilin National High-tech Zone Information Industrial Park,
Chaoyang Road, Qixing District, Guilin, Guangxi, 541004, P.R. China

EC	REP	Xpedent Austria GmbH Millenium Park 9, 6890 Lustenau, Austria
-----------	------------	---

POLO | DENT[®]

EQUIPMENT

Product #	Description	Product #	Description
	Scaling Tips:	52 04 40	(Sirona) no. PD-1L
52 03 00	(EMS) no. PD-A (G4)	52 04 45	(Sirona) no. PD-2L
52 03 10	(EMS) no. PD-P (G3)	52 04 50	(Sirona) no. PD-3L
52 03 20	(EMS) no. PD-PS (P3)	52 04 55	(Sirona) no. PD-4L
52 03 40	(KaVo) no. PD-5 (GK1#5)	52 04 60	(Sirona) no. PD-PE1
52 03 45	(KaVo) no. PD-6 (GK2#6)	52 04 65	(Sirona) no. PD-PE2
52 03 50	(KaVo) no. PD-7 (GK3#7)	52 04 70	(Sirona) no. PD-PE3
52 03 55	(KaVo) no. PD-8 (GK4#8)		
52 03 58	(KaVo Quick) no.PD-5A		Torque Wrench for Scaling Tips:
52 03 59	(KaVo Quick) no.PD-6A	52 04 80	(EMS)
52 03 60	(KaVo Quick) no.PD-7A	52 04 85	(KaVo)
52 03 61	(KaVo Quick) no.PD-8A	52 04 90	(Acteon)
52 03 65	(KaVo) no.PD-9 (GC2)	52 04 95	(Sirona)
52 03 66	(KaVo) no.PD-10 (GC1)	52 05 00	(Endo)
52 03 67	(KaVo) no.PD-11 (PC1)		
52 03 69	(KaVo Piezo) no.PD-201		Implant Cleaning Tips Set:
52 03 70	(KaVo Piezo) no.PD-202	52 03 25	(EMS) no.PD-eICP SET
52 03 71	(KaVo Piezo) no.PD-203	52 03 74	(KaVo) no.PD-kICP SET
52 03 75	(Acteon) no. PD-1	52 04 20	(Acteon) no.PD-slCP SET
52 03 80	(Acteon) no. PD-2	52 04 75	(Sirona) no.PD-srICP SET
52 03 85	(Acteon) no. PD-1S		
52 03 90	(Acteon) no. PD-10P		
52 03 95	(Acteon) no. PD-10X		
52 04 00	(Acteon) no. PD-10Z		
52 04 05	(Acteon) no. PD-sH3		

Scaling Tips

ENGLISH

POLO | DENT |[®]
EQUIPMENT

Indications for Use:

Scaling tips work with ultrasonic scaler machines. The ultrasonic energy is created by the scaler machine and converted into vibrations by the handpiece; these vibrations are then focused and applied to the area to be treated by the scaling tip. Ultrasonic scaler machines and tips are only to be used by qualified dental professionals, either by, or on the order of a dentist. Tips are divided into the following categories: General (scaling-prophylaxis), Periodontal, Endodontic, Prosthesis, Implantation, Sinus Lifting, Extraction, Cutting.

General (Scaling–prophylaxis)

These are used for the removal of supragingival stains and deposits. The power setting used should be the lowest effective power according to the instructions of the scaler machine.

Periodontal

These are used for subgingival work. The power settings will depend upon the tip used and the operation to be carried out; in general though, thinner tips require lower power.

Endodontic

These are used for endodontic treatments and are all longer, thinner tips. Endo- use will either have a special power setting or require a low power to be selected.

Prosthesis

The removal of crowns or similar prosthetics may require a higher power setting.

Implantation

Used for preparation and optimization of implant site.

Sinus Lifting

Used for preparation, separation and lifting of the sinus membrane.

Extraction

Used for root osteoplasty.

Cutting

Used for root preparation, chip harvesting, osteotomy and osteoplasty.

General Instructions:

1. Contraindications: Follow scaler machine instructions for contraindications.
2. Scaler machine's instructions: Please read and follow your scaler machine's instructions.
3. Cleaning and sterilization: Scaling tips must be sterilized before first use and after each use (they are not delivered sterile).

Cleaning

The removal of organic deposits is very important and will achieve best results by putting them immediately after use into 70% Isopropyl alcohol (Isopropanol) or a low-alcohol cleaning solution. After that use an ultrasonic bath with a low-foaming detergent. Alternatively, you can clean them in warm water with a low-foaming detergent and an autoclavable brush.

Attention

Do not use chemical disinfectants as these can damage the scaling tips.

Scaling Tips

ENGLISH

POLO | DENT |[®]
EQUIPMENT

After cleaning

Rinse the cleaned tips thoroughly in warm water to remove all chemicals and let them dry completely. PoloDent scaling tips are made out of stainless steel of the highest quality. Therefore, the use of an antirust agent before steam sterilization is acceptable, although it is not required. Pack them in suitable packaging, bags, trays or cassettes. Check them with spore test sets. Biological monitoring devices should be used at least once a week. In addition, every time if new packaging, new staff, new procedures, new equipment or repaired equipment is used, a biological monitoring should be added to check if these changes are still effective.

Attention

Do not use closed containers during sterilization with steam or ethylene oxide. The cleaning and disinfecting device must be approved by its manufacturer for cleaning and disinfection be released by dental instruments and EN ISO 15883-1.

Sterilize with ethylene oxide

Put the dried packaged materials in the unit and work according to the manufacturer's instructions.

Sterilize with steam

Place the dried wrapped tip in an autoclave. Recommended setting: 134°C for 18 minutes.

Regarding sterilization

Ethylene oxide is the preferred method for sterilization. Steam autoclaving is the preferred method of dentistry. Dry heat should NOT be used for materials which can be sterilized by steam. Chemical agents are not recommended since their effectiveness cannot be monitored biologically, their cycle times are longer (6-10 hours), and all rinsing and handling must be aseptic. Dry Heat and chemical steam sterilization are considered ineffective methods, with the risk of damaging materials. DO NOT use chemiclaves or drying ovens for sterilizing PoloDent scaling tips.

The scaling tip has no sterilization cycle limit, but after the tip is worn exceeding 2 mm, the scaling tip needs to be replaced.

4. Tip selection: Standard tips are best for medium to heavy deposits. Universal tips are designed for use in pockets less than 4mm.
5. Tip attachment: Tips should be loaded onto the handpiece using the appropriate Torque Wrench.
6. Power setting: When the scaler tip is used with the handpiece for the first time, set the scaler machine to the minimum power setting, and then adjust to the correct power, according to the machine's instructions. The use of excessive power may cause patient discomfort and may cause micro fractures in the tip.
7. Cooling: In order to cool the scaling tip in use, the power and then the water flow need to be set, follow machine instructions for this;
8. Tip modification: DO NOT modify, sharpen or bend the tip in any way; this will result in incorrect ultrasonics and, will damage (or possibly break) the tip and the scaler machine.
9. Damaged tips: If tip is bent out of shape, if a tip has been dropped, or if a change of power is experienced during use, discard and replace.
10. Storage environment: Scaling tips should be autoclaved after use. Scaler tips should be stored in a dry environment with humidity lower than 70%.
11. Tip Life: Scaling tips effectiveness is reduced by wear; once the tip has worn by 2 mm, it will no longer be effective. As a guide, basic scaling tips should last for about 100 uses, with Endo tips needing much more frequent change; thinner tips have shorter life. Some of the specialist tips, such as those with diamond or Titanium Nitride coats will be replaced after just 3 or 4 uses.

Any serious incident related to this medical device must be reported to the manufacturer and the competent authorities.

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